

Datalogging experiment (4)



**ASSIGN
BUSTER**

Datalogging Experiment (4) Acid-base Titration using Method of Double

Indicators Student Handout Purposes To determine the composition of the following mixture by double indicator method: 1. NaOH(aq) and Na₂CO₃(aq)

2. NaHCO₃(aq) and Na₂CO₃(aq) Introduction Consider a mixture of

NaOH(aq) and Na₂CO₃(aq). Reaction between HCl(aq) and Na₂CO₃(aq)

takes place in two stages: $\text{HCl(aq)} + \text{Na}_2\text{CO}_3\text{(aq)} \rightarrow \text{NaHCO}_3\text{(aq)} +$

$\text{H}_2\text{O(l)}$ (1) $\text{HCl(aq)} + \text{NaHCO}_3\text{(aq)} \rightarrow \text{NaCl(aq)} + \text{CO}_2\text{(g)} +$

$\text{H}_2\text{O(l)}$ (2) While that between HCl(aq) and NaOH(aq) completes in

only one step: $\text{HCl(aq)} + \text{NaOH(aq)} \rightarrow \text{NaCl(aq)} + \text{H}_2\text{O(l)}$

..... (3) Solution mixture of reaction (1) at the equivalence point is

alkaline, that of reaction (2) is acidic and that of reaction (3) is neutral. Thus

the whole titration should have three breaks in the pH curve, corresponding

to the above three stages. Reactions (1) and (3) can be indicated by

phenolphthalein and that of reaction (2) can be indicated by methyl orange.

Stoichiometry confines each of the above pH reactions to react according to

a mole ratio of 1 : 1. This means, say from equation (2), the number of mole

of HCl(aq) determined from the methyl orange titration is equal to the

number of mole of NaHCO₃(aq). Likewise, total number of moles of

NaOH(aq) and Na₂CO₃(aq) in the solution mixture can be calculated

according to the volumes of HCl(aq) added at the end point Vol. of HCl

indicated by the colour change of the phenolphthalein indicator.

Alternatively, the Fig. 1: Titration curve for a mixture of three break points

(see Fig. 1) also indicate NaOH(aq) and Na₂CO₃(aq) with HCl(aq) the volume

of HCl(aq) required for each reaction. 49 Datalogging Experiment (4) For the

titration of a mixture of NaHCO₃(aq) pH and Na₂CO₃(aq) with HCl(aq), only

two break points are expected (see Fig. 2). Volume of HCl(aq) added for each

break point can be easily obtained by observing either the colour change at the end point or the shape of the titration curve. Vol of HCl Fig. 2: Titration

curve for a mixture of $\text{NaHCO}_3(\text{aq})$ and $\text{Na}_2\text{CO}_3(\text{aq})$ with $\text{HCl}(\text{aq})$ Safety

Avoid skin contact with chemicals. Any acid or alkali spilt should be

thoroughly washed with tap water. Materials and Apparatus Phenolphthalein

indicator, methyl orange indicator FLAMMABLE EYE PROTECTION MUST BE

WORN Standard 0.15 M and 0.2 M $\text{HCl}(\text{aq})$, a mixture of $\text{NaOH}(\text{aq})$ and

$\text{Na}_2\text{CO}_3(\text{aq})$, a mixture of $\text{NaHCO}_3(\text{aq})$ and $\text{Na}_2\text{CO}_3(\text{aq})$, deionised water, a

datalogger with pH sensor, computer, magnetic stirrer, small beaker, burette

and pipette, stand and burette clamp. Experimental Procedures Part A:

Titration of a mixture of $\text{NaOH}(\text{aq})$ and $\text{Na}_2\text{CO}_3(\text{aq})$ with 0.15 M $\text{HCl}(\text{aq})$

using phenolphthalein indicator followed by methyl orange indicator 1. 2. Set

up the interface box and connect it to the computer. Arrange the setup for

pH determination. The pH sensor should be calibrated before use. Pipette 25

cm^3 of the solution mixture into a small beaker and add 2 drops of

phenolphthalein indicator. Place a stirrer bar into the alkaline solution and

rest the beaker on a magnetic stirrer which is covered by a white tile. Switch

on the magnetic stirrer. Lower the pH electrode into the alkaline solution,

ensuring that the glass bulb is completely immersed while the stirrer bar is

spinning smoothly (see Fig. 3). Fig. 3: Setup for pH datalogging 50

Datalogging Experiment (4) 3. Start the datalogging software and select the

pH sensor with graph display screen format. Set the pH value limits from 0 to

13. The logging mode is set to manual operation with readings taken for

each cm^3 addition of the titrant. Fill the burette with 0.15 M $\text{HCl}(\text{aq})$. Turn

the stopcock of the burette open and start recording at the same time. When

approaching the first end point, as judged by the rate of disappearance of

the red colour of the solution in the beaker, adjust interval additions each to 0.50 cm³. Allow time for the evolution of CO₂ gas bubbles as their presence round the glass bulb of the electrode may interfere with pH measurements. Observe the colour change of the reaction mixture and the pH values displayed on the computer screen carefully when the titration passes the first end point. Add 3 drops of methyl orange indicator when the reaction mixture becomes colourless. Look for the second end point and continue titration as in step (4) until a total volume of 50 cm³ of the titrant has been added. Save the data file. Carefully empty the contents of the beaker, pay special attention to retain the small stirrer bar. Clean the pH electrode with deionised water.

4. 5. 6. 7. 8. Part B 9. As described above, titrate a mixture of NaHCO₃(aq) and Na₂CO₃(aq) with 0.20 M HCl(aq) using phenolphthalein indicator followed by methyl orange indicator. Results Table A

Phenolphthalein indicator	Final burette reading/cm ³	Initial burette reading/cm ³	Volume of 0.15 M HCl used/cm ³
51			

Table B Phenolphthalein indicator

Phenolphthalein indicator	Final burette reading/cm ³	Initial burette reading/cm ³	Volume of 0.20 M HCl used/cm ³
51			

Datalogging Experiment (4) Treatment of Data Part A 1. From the methyl orange end point, calculate the number of moles of 0.15 M HCl(aq) added and hence the number of moles of Na₂CO₃(aq) in 25 cm³ of the alkaline solution mixture. From the phenolphthalein end point, calculate the number of moles of 0.15 M HCl(aq) added and hence the total number of moles of NaOH(aq) and Na₂CO₃(aq) in 25 cm³ of the solution mixture. Calculate the number of mole of NaOH(aq) in 25 cm³ of the solution mixture. Calculate the mass of Na₂CO₃ and NaOH in 1 dm³ of the solution mixture respectively.

2. 3. 4. Part B 1. From the methyl orange end point, calculate the number of

moles of 0.2 M HCl(aq) added and hence the number of moles of Na₂CO₃(aq) in 25 cm³ of the solution mixture. From the phenolphthalein end point, calculate the number of moles of 0.2 M HCl(aq) added and hence the number of moles of NaHCO₃(aq) in 25 cm³ of the solution mixture.

Calculate the mass of Na₂CO₃ and that of NaHCO₃ in 1 dm³ of the solution mixture.

2. 3. Discussion Questions

1. 2. Suggest other indicators that can be used in place of methyl orange and phenolphthalein. Explain. Can the same method be applied to determine the concentrations of Na₃PO₄(aq) and NaH₂PO₄(aq) in a solution mixture of the two salts? What factors should be considered?

52